

Revaclear

DESIGNED FOR:
HFHD [High flux]

OTHER APPLICABLE THERAPIES:
CONVECTIVE [HDF]

MEMBRANE:
PORACTON [PAES/PVP, BPA-free]

OPTIMIZING HIGH FLUX PERFORMANCE

The **Revaclear** dialyzer series is a range of high efficiency high-flux dialyzers designed to enhance safety and biocompatibility for your patients, while optimizing clearance with a smaller surface area.¹

OPTIMIZING PERFORMANCE FOR ALL YOUR PATIENTS²

The **Revaclear** dialyzers are designed to optimize the performance of high-flux treatments.

- The **Poracton** membrane provides effective permeability with minimal resistance to diffusion^{3,4,5}
- Three surface area options are available to meet individual patient needs
- Study in HD showed **Revaclear 400** to remove small solutes and β_2 -microglobulin to a similar extent as a 22% larger surface area dialyzer²

WITH SAFETY AND BIOCOMPATIBILITY IN MIND

The relative compact surface area of the **Revaclear** dialyzers may help manage some patient risks.

- Reduces exposure of blood, potentially reducing clotting and micro inflammation⁶
- Produces less biohazardous waste and reduces saline need, compared to dialyzers of the same performance^{7,8}



Revaclear Specifications

| MATERIALS | REVACLEAR 300 | REVACLEAR 400 | REVACLEAR 500 |
|-----------------|--|---------------|---------------|
| Membrane | Poracton Polyarylethersulfone and Polyvinylpyrrolidone blend BPA-free | | |
| Potting | Polyurethane (PUR) | | |
| Housing | Polycarbonate (PC) | | |
| Gaskets | Silicone rubber (SIR) | | |
| Protection caps | Polypropylene (PP) | | |
| Sterilization | Steam (inside-out) | | |
| Sterile barrier | Tyvek | | |

| SPECIFICATIONS | | | |
|---|------------------|---------|---------|
| UF-Coefficient (mL/h*mmHg)* | 48 | 54 | 65 |
| KoA urea* | 1186 | 1439 | 1578 |
| Blood Compartment volume (mL) | 74 | 93 | 106 |
| Minimum recommended priming volume (mL) | 300 | | |
| Maximum TMP (mmHg) | 600 | | |
| Recommended Q _B (mL/min) | 200-500 | 200-600 | 250-600 |
| Storage conditions | <30°C (or <86°F) | | |
| Units per box | 24 | | |
| Gross/net weight (g) | 215/160 | 225/170 | 250/190 |

| MEMBRANE | | | |
|---|-----|-----|-----|
| Effective Membrane Area (m ²) | 1.4 | 1.8 | 2.1 |
| Fiber inner diameter (µm) | 190 | | |
| Fiber wall thickness (µm) | 35 | | |

| SIEVING COEFFICIENTS* | |
|--|--------|
| Vitamin B12 (1,4 kDa) | 1.0 |
| Inulin (5,2 kDa) | 1.0 |
| β ₂ -microglobulin (11,8 kDa) | 0.95 |
| Myoglobin (17 kDa) | 0.68 |
| Albumin (66,4 kDa) | 0.0027 |

* According to EN 1283/ISO 8637:
 - UF-Coefficient: measured with bovine blood, Hct 32%, Pct 60g/L, 37°C
 - KoA urea: calculated at Q_B=300 mL/min, Q_D=500mL/min, UF=0 mL/min
 - Sieving coefficients: measured with human plasma, Q_B=300 mL/min, UF=60 mL/min
 - Clearances In-Vitro: measured at UF=0 mL/min, ±10%

| CLEARANCES IN VITRO (mL/min)* | REVACLEAR 300 | REVACLEAR 400 | REVACLEAR 500 |
|---|---------------|---------------|---------------|
| Urea (60 Da) (Q_B-Q_D, mL/min) | | | |
| 200-250**/500 | 196 | 198 | 244 |
| 300/500 | 272 | 281 | 284 |
| 400/500 | 323 | 338 | 345 |
| 400/800 | 355 | 369 | 375 |
| 500/800 | 408 | 430 | 439 |
| Phosphate (95 Da) | | | |
| 200-250**/500 | 191 | 195 | 238 |
| 300/500 | 256 | 267 | 272 |
| 400/500 | 298 | 315 | 323 |
| 400/800 | 330 | 348 | 355 |
| 500/800 | 373 | 398 | 409 |
| Creatinine (113 Da) | | | |
| 200-250**/500 | 185 | 191 | 230 |
| 300/500 | 242 | 255 | 261 |
| 400/500 | 278 | 297 | 306 |
| 400/800 | 309 | 330 | 338 |
| 500/800 | 345 | 373 | 384 |
| Vitamin B12 (1.4 kDa) | | | |
| 200-250**/500 | 146 | 158 | 183 |
| 300/500 | 174 | 191 | 200 |
| 400/500 | 191 | 213 | 223 |
| 400/800 | 212 | 236 | 247 |
| 500/800 | 228 | 256 | 269 |

** REVACLEAR 500

1. Baxter. *REVACLEAR White Paper*. USMP/MG3/140052, May 2013.
2. Mauric A, et al. *Poster SP401*, presented at 50th ERA-EDTA congress. Istanbul (Turkey), 2013.
3. Ronco C, et al. *Evolution of synthetic membranes for blood purification: the case of the Polyflux family*. *Nephrol Dial Transplant* 2003;18(Suppl 7):vii10-20.
4. Ward R, et al. *Abstract SA-P0510*, presented at the 40th ASN congress. San Francisco (USA), 2007.
5. Bhimani JP, et al. *Effect of increasing dialysate flow rate on diffusive mass transfer of urea, phosphate and beta2-microglobulin during clinical haemodialysis*. *Nephrol Dial Transplant* 2010; 25:3990-3995.
6. Yao Q, et al. *Inflammation as a cause of malnutrition, atherosclerotic cardiovascular disease, and poor outcome in hemodialysis patients*. *Hemodial Int* 2004; 8:118-129.
7. Baxter. Data on file. *Biohazardous waste cost calculation*, 2015.
8. Baxter. *REVACLEAR dialyzer priming guide*. 306150152_C, 2009.

The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (Full quality assurance system of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/47/EC)

For safe and proper use of the device, please refer to the Instructions for Use

