



# The role of low-flow ECCO<sub>2</sub>R in supporting LPV strategies

**Baxter**





# The role of low-flow ECCO<sub>2</sub>R in supporting LPV strategies



How may mechanical ventilation adversely affect patient outcomes?



How may LPV reduce risk of VILI?



How may ECCO<sub>2</sub>R facilitate the use of LPV?

**Abbreviations:** ALI, acute lung injury; ARDS, acute respiratory distress syndrome; CI, confidence interval; Crs, respiratory system compliance; ECCO<sub>2</sub>R, extracorporeal carbon dioxide removal; FiO<sub>2</sub>, fraction of inspired oxygen; ICU, intensive care unit; LPV, lung protective ventilation; MV, mechanical ventilation; OR, odds ratio; PaCO<sub>2</sub>, arterial carbon dioxide partial pressure; PaO<sub>2</sub>, arterial oxygen partial pressure; PBW, predicted body weight; PEEP, positive end-expiratory pressure; P<sub>plat</sub>, plateau pressure; RR, relative risk; V<sub>A</sub>/Q, alveolar ventilation to perfusion ratio; VILI, ventilator-induced lung injury; V<sub>T</sub>, tidal volume; VV-ECMO, veno-venous extracorporeal membrane oxygenation; ΔP, driving pressure.



# How may mechanical ventilation adversely affect patient outcomes?


## VILI is a potential complication of mechanical ventilation

- The goal of mechanical ventilation is to provide acceptable oxygenation and CO<sub>2</sub> removal while minimizing VILI<sup>1</sup>
- Complications of mechanical ventilation include volutrauma, barotrauma, atelectrauma, and biotrauma<sup>2,3</sup>

### Complications of mechanical ventilation

<b>Volutrauma (biophysical injury)</b>
Over-distension of alveoli resulting from increased $V_T$
<b>Barotrauma (biophysical injury)</b>
Alveolar rupture and air leaks resulting from high pressure
<b>Atelectrauma (biophysical injury)</b>
Damage caused by repetitive opening and closing of collapsed lung parts
<b>Biotrauma (biochemical injury)</b>
Release of biological mediators and translocation into the circulation

### Mechanisms of VILI



Adapted from Gattinoni L, Protti A. *CMAJ* 2008;178:1174–6, with permission from Access Copyright.

- ARDS accounts for approximately 25% of patients requiring mechanical ventilation<sup>4</sup>
- As most patients with ARDS require invasive mechanical ventilation,<sup>1</sup> they are at risk of VILI





## ARDS

- An acute inflammatory lung injury that leads to increased pulmonary vascular permeability, increased lung weight and loss of aerated lung tissue, resulting in hypoxemia and bilateral radiographic opacities<sup>5</sup>
- A serious condition that is common but often under-recognized in the ICU,<sup>4</sup> which may limit implementation of effective management
  - In the LUNG SAFE study, 35.8% of all cases were not recognized by physicians<sup>4</sup>
  - Less severe ARDS was more likely to be unrecognized (48.7% of mild cases vs 21.5% of severe cases)<sup>4</sup>



### High incidence:<sup>4</sup>

- 1 in 10 of all ICU admissions\*
- 1 in 4 of all patients requiring mechanical ventilation\*



### High mortality rate:<sup>4</sup>

- 40.0% hospital mortality\*
- 35.3% ICU mortality\*



### Risk of multiple organ failure:

- 68.9–70.0% risk of failure of  $\geq 2$  organs<sup>6</sup>
- 30.0–31.1% risk of failure of  $\geq 3$  organs<sup>6</sup>
- High rates of renal failure (41–49%) and liver failure (13–34%)<sup>7,8</sup>

\*As reported in LUNG SAFE; a large, international, prospective, cohort study (n = 2377).<sup>4</sup>

## Evidence from animal studies suggests that VILI may contribute to the development of multiple organ failure

- A proposed mechanism for a relationship with multiple organ failure is based on the systemic release of inflammatory mediators resulting from VILI (biotrauma)<sup>3,9</sup>
- Multiple organ failure has been associated with increased risk of mortality in patients with ALI or ARDS; in one study, multiple organ failure was the cause of death in 16.7% of patients with ARDS<sup>6</sup>

## Key points



- Mechanical ventilation is the cornerstone of treatment for patients with impaired lung function<sup>1</sup>
- However, VILI may complicate the management of mechanically ventilated patients,<sup>1,2</sup> particularly those with ARDS
- Evidence from animal studies suggests that VILI can contribute to poor outcomes, including multiple organ failure<sup>3,9</sup>





# How may LPV reduce risk of VILI?

## Lung protective ventilation strategies modify ventilation parameters to reduce the risk of VILI

<b>Conventional MV</b>	<ul style="list-style-type: none"> <li><math>V_T</math> of 10–15 mL/kg PBW has been traditionally used to normalize <math>\text{PaCO}_2</math>, <math>\text{PaO}_2</math> and pH<sup>9</sup></li> <li>May exacerbate or perpetuate lung injury<sup>9</sup></li> </ul>
<b>Concept of LPV</b>	<ul style="list-style-type: none"> <li>Utilizes lower <math>V_T</math> (~6 mL/kg PBW) than conventional MV<sup>9</sup></li> <li>Other components may include lower <math>P_{\text{plat}}</math>, higher PEEP, and lower <math>\Delta P</math><sup>8,10,11</sup></li> <li>Elevated <math>\text{PaCO}_2</math> is either accepted (permissive hypercapnia) or may require measures to reduce <math>\text{CO}_2</math> levels<sup>12</sup></li> </ul>
<b>Concept of ultra LPV</b>	<ul style="list-style-type: none"> <li>Utilizes even lower <math>V_T</math> (<math>\leq 3</math> mL/kg PBW) compared with LPV<sup>13</sup></li> <li>Other components may include lower <math>P_{\text{plat}}</math>, higher PEEP, and lower <math>\Delta P</math><sup>14</sup></li> <li>Greater elevations in <math>\text{PaCO}_2</math> can occur compared with LPV and extracorporeal lung support is needed to reduce <math>\text{CO}_2</math> levels<sup>13</sup></li> </ul>



### Driving pressure

$\Delta P$  is defined as  $V_T$  normalized to Crs ( $V_T/\text{Crs}$ ) or  $P_{\text{plat}}$  minus PEEP

$$\Delta P = P_{\text{plat}} - \text{PEEP}$$

## LPV reduces duration of mechanical ventilation and mortality compared with conventional ventilation


The landmark, randomized, controlled ARDS Network (ARDSNet) study compared conventional ventilation with LPV in 861 patients with ARDS<sup>9</sup>

- Conventional ventilation:** initial  $V_T$  12 mL/kg PBW and  $P_{\text{plat}} \leq 50$  cmH<sub>2</sub>O
- LPV:** initial  $V_T$  6 mL/kg PBW and  $P_{\text{plat}} \leq 30$  cmH<sub>2</sub>O

The LPV strategy was associated with:<sup>9</sup>

- Greater probability of survival over 180 days (69.0% vs 60.2%;  $p = 0.007$ )
- Greater number of ventilator-free days during the first 28 days (12 vs 10 days;  $p = 0.007$ )
- Greater number of non-pulmonary organ/system failure-free days during the first 28 days (15 vs 12 days;  $p = 0.006$ )

## LPV with the ARDSNet protocol reduces mortality compared with conventional mechanical ventilation<sup>9</sup>




Mortality rate was reduced by 8.8 percentage points with the LPV protocol vs a conventional ventilation protocol

(31.0% vs 39.8%;  $p = 0.007$ )

**Randomized controlled trial (ARDSNet):** Clinical outcomes of patients with ALI or ARDS managed with a LPV protocol (initial  $V_T$  6 mL/kg PBW and  $P_{plat}$  maintained between 25–30 cmH<sub>2</sub>O) were compared with those managed with a conventional ventilation protocol (initial  $V_T$  12 mL/kg PBW and  $P_{plat} \leq 50$  cmH<sub>2</sub>O). Adapted with permission from ARDS Network. *N Engl J Med* 2000;342:1301–8. Copyright © (2000) Massachusetts Medical Society.

## Reducing driving pressure as part of an LPV strategy may improve survival<sup>11</sup>



•  $\Delta P$  was strongly associated with mortality

(RR 1.41; 95% CI 1.31–1.51,  $p < 0.001$ )

• This was true even in patients receiving 'protective'  $P_{plat}$  and  $V_T$

(RR 1.36; 95% CI 1.17–1.58,  $p < 0.001$ )

• Reductions in  $V_T$  or increases in PEEP were associated with better survival only if associated with reductions in  $\Delta P$

Median  $V_T$  (10<sup>th</sup>–90<sup>th</sup> percentile) mL/kg PBW

6.0 (5.9–7.5)

6.1 (5.8–9.2)

8.0 (5.7–12.1)

**Randomized controlled trials:** The relationship between different ventilation parameters and mortality was explored in a mediation analysis of data from 9 randomized controlled trials in patients with ARDS ( $n = 3562$ ). Data shown as the increase in RR of hospital mortality as a function of  $\Delta P$  after multivariate adjustment (95% CIs represented as grey shaded area). Adapted with permission from Amato MB, et al. *N Engl J Med* 2015;372:747–55. Copyright © (2015) Massachusetts Medical Society.

## Key points

- LPV strategies modify ventilation parameters that have been shown to increase the risk of VILI<sup>9</sup>
- The ARDSNet LPV protocol, based on reduced  $V_T$  and  $P_{plat}$ , improves patient outcomes, including mortality<sup>9,11</sup> and duration of mechanical ventilation<sup>9</sup>
- More recent evidence shows that lower  $\Delta P$  is associated with a reduced mortality risk in patients with ARDS<sup>11</sup>






# How may low-flow ECCO<sub>2</sub>R facilitate the use of LPV?

## Despite guidelines supporting the use of LPV,<sup>15–18</sup> $V_T$ often exceeds 6 mL/kg PBW in clinical practice

- In a cross-sectional survey of 200 German ICUs ( $n = 152$ ), only 2.6% of patients received low  $V_T$  ventilation despite the fact that perceived adherence by ICU directors was 79.9%<sup>19</sup>
- In the LUNG SAFE study ( $n = 2255$ ), more than one-third of patients were mechanically ventilated with  $V_T > 8$  mL/kg PBW<sup>4</sup>

### Adherence to LPV is poor in clinical practice<sup>4</sup>




- $V_T \leq 6$  mL/kg PBW was used in only about 20% of patients with ARDS in the LUNG SAFE study
- Differences in mean  $V_T$  between patients with mild and severe ARDS were clinically modest

**Prospective cohort study:** Ventilatory management of patients with ARDS (Berlin definition;  $n = 2255$ ) was assessed as a secondary endpoint in the LUNG SAFE study. Adapted with permission from Bellani G, et al. *JAMA* 2016;315:788–800. Copyright © (2016) American Medical Association. All rights reserved.

## Major barriers to LPV adherence include concerns about hypercapnia and respiratory acidosis induced by $V_T$ reduction<sup>20–22</sup>

- Hypercapnia is often regarded as an acceptable side effect, however, physiological effects may include pulmonary vasoconstriction, increased intracranial pressure, and decreased renal blood flow, among others<sup>12</sup>
- Recent evidence indicates that hypercapnia is associated with an increased risk of ICU mortality<sup>23</sup>

## Severe hypercapnia is associated with an increased risk of ICU mortality<sup>23</sup>



PaCO<sub>2</sub> ≥ 50 mmHg was independently associated with higher ICU mortality (OR 2.40; 95% CI 1.67–3.46, p < 0.001)

**Prospective observational studies:** The relationship between hypercapnia and ICU mortality was assessed in a secondary analysis of 3 studies that included data from 1899 patients with moderate-to-severe ARDS. Results from a logistic regression model with adjustment for baseline variables are shown. Adapted with permission of Springer from Nin N, et al. *Intensive Care Med* 2017;43:200–8.

## ECCO<sub>2</sub>R enhances CO<sub>2</sub> removal in patients receiving LPV

- Evidence shows that ECCO<sub>2</sub>R systems significantly reduce PaCO<sub>2</sub> levels, and may, therefore, facilitate LPV by allowing for a reduction in V<sub>T</sub><sup>24</sup>
- More specifically, low-flow ECCO<sub>2</sub>R devices using flow rates as low as 0.5 L/min should theoretically be sufficient to eliminate all CO<sub>2</sub> produced by the body<sup>24</sup>




1 L of blood with a PaCO<sub>2</sub> of 5 kPa contains around 500 mL of CO<sub>2</sub>, or on average, two times more CO<sub>2</sub> than the body produces per minute

## Low-flow ECCO<sub>2</sub>R enables use of LPV by reducing PaCO<sub>2</sub> levels and normalizing arterial pH

- The ability of low-flow ECCO<sub>2</sub>R to facilitate use of LPV has been demonstrated in a prospective study of patients with ARDS, 10 of whom had P<sub>plat</sub> within the range of 28–30 cmH<sub>2</sub>O while being treated with the ARDSNet protocol<sup>25</sup>



## Low-flow ECCO<sub>2</sub>R allows for maintenance of $V_T$ and $P_{plat}$ in line with LPV strategies<sup>25</sup>




**Prospective study:** Data shown as individual and average (horizontal bars) values of  $V_T$  and  $P_{plat}$  during LPV with the ARDSNet protocol, after lowering  $V_T$  (BL), and at 1–1.5, 24, 48, and 72 hours after initiation of ECCO<sub>2</sub>R (n = 10 ICU patients with ARDS [American-European Consensus Conference definition]). \*p < 0.001 vs ARDSNet ventilation. BL, baseline. Adapted with permission from Terragni PP, et al. *Anesthesiology* 2009;111:826–35.



- Patients received low-flow ECCO<sub>2</sub>R following a reduction in  $V_T$  to < 6 mL/kg PBW
- Low-flow ECCO<sub>2</sub>R allowed for maintenance of  $V_T$  < 6 mL/kg/PBW
- Low-flow ECCO<sub>2</sub>R allowed for maintenance of  $P_{plat}$  between 25–28 cmH<sub>2</sub>O

## Low-flow ECCO<sub>2</sub>R reduces PaCO<sub>2</sub>, thereby normalizing arterial pH<sup>25</sup>




**Prospective study:** Data shown as individual and average (horizontal bars) values of PaCO<sub>2</sub> and arterial pH during LPV with the ARDSNet protocol, after lowering  $V_T$  (BL), and at 1–1.5, 24, 48, and 72 hours after initiation of ECCO<sub>2</sub>R (n = 10 ICU patients with ARDS [American-European Consensus Conference definition]). \*p < 0.001 vs ARDSNet ventilation. Adapted with permission from Terragni PP, et al. *Anesthesiology* 2009;111:826–35.



- Compared with baseline, PaCO<sub>2</sub> was significantly reduced and arterial pH was significantly increased after 1–1.5 hours of initiating low-flow ECCO<sub>2</sub>R (p < 0.001 for both)
- After 72 hours of low-flow ECCO<sub>2</sub>R, the reduction in PaCO<sub>2</sub> was sufficient to normalize pH


## Ultra-protective ventilation facilitated by low-flow ECCO<sub>2</sub>R may reduce driving pressure<sup>14</sup>



- ΔP – an independent risk factor for hospital mortality<sup>11</sup> – was significantly reduced during the first two days of treatment vs baseline (13.9 vs 11.6 cmH<sub>2</sub>O;  $p < 0.05$ )
- PaCO<sub>2</sub> and pH were corrected to within 10% of baseline values

**Prospective study:** The feasibility of very low  $V_T$  ventilation (4 mL/kg PBW) combined with low-flow ECCO<sub>2</sub>R was evaluated in patients with moderate ARDS (Berlin definition;  $n = 15$ ). Adapted from Fanelli V, et al. *Crit Care* 2016;20:36.

## ECCO<sub>2</sub>R may have a role in the management of patients with moderate to severe ARDS – an example algorithm<sup>26</sup>



**Assuming:** Optimum MV settings and application of adjunctive therapies

Assuming: Optimum MV settings and application of adjunctive therapies		
<b>Rare circumstances:</b> ECCO <sub>2</sub> R may be considered to facilitate ultra LPV	<b>Consider ECCO<sub>2</sub>R:</b> <ul style="list-style-type: none"> <li>• pH &lt; 7.25 for &gt; 2 h</li> <li>• <math>P_{\text{plat}} &gt; 28</math> cmH<sub>2</sub>O for &gt; 2 h</li> <li>• Further reduction in risk of VILI (ultra LPV)</li> </ul>	<b>Consider VV-ECMO to facilitate gas exchange during complete lung rest*</b>
<b>Mild ARDS</b>	<b>Moderate ARDS</b>	<b>Severe ARDS</b>

Adapted from Del Sorbo L, et al. *Lancet Respir Med* 2014;2:154–64, with permission from Elsevier.

**Note:** This example algorithm is based on author opinion and is not a recognized guideline. Severity based on the Berlin definition.<sup>5</sup>

\*PaO<sub>2</sub>/FiO<sub>2</sub> < 50 with FiO<sub>2</sub> > 0.8 for > 3 h, or PaO<sub>2</sub>/FiO<sub>2</sub> < 80 with FiO<sub>2</sub> > 0.8 for > 6 h.

### Key points

- Despite evidence and guidelines in support of LPV, concerns about hypercapnia and respiratory acidosis resulting from LPV may limit its use<sup>20–22</sup>
- Recent evidence shows that severe hypercapnia is independently associated with increased ICU mortality<sup>23</sup>
- ECCO<sub>2</sub>R – even at low blood flow rates – enables use of LPV and ultra LPV by reducing PaCO<sub>2</sub> and thereby normalizing pH<sup>25</sup>





# The role of low-flow ECCO<sub>2</sub>R in supporting LPV strategies



## How may mechanical ventilation adversely affect patient outcomes?

- Mechanical ventilation is the cornerstone of treatment for patients with impaired lung function<sup>1</sup>
- However, VILI may complicate the management of mechanically ventilated patients,<sup>1,2</sup> particularly those with ARDS
- Evidence from animal studies suggests that VILI can contribute to poor outcomes, including multiple organ failure<sup>3,9</sup>



## How may LPV reduce risk of VILI?


- LPV strategies modify ventilation parameters that have been shown to increase the risk of VILI<sup>9</sup>
- The ARDSNet LPV protocol, based on reduced  $V_T$  and  $P_{plat}$ , improves patient outcomes, including mortality<sup>9,11</sup> and duration of mechanical ventilation<sup>9</sup>
- More recent evidence shows that lower  $\Delta P$  is associated with a reduced mortality risk in patients with ARDS<sup>11</sup>



## How may low-flow ECCO<sub>2</sub>R facilitate the use of LPV?

- Despite evidence and guidelines in support of LPV, concerns about hypercapnia and respiratory acidosis resulting from LPV may limit its use<sup>20–22</sup>
- Recent evidence shows that severe hypercapnia is independently associated with increased ICU mortality<sup>23</sup>
- ECCO<sub>2</sub>R – even at low blood flow rates – enables use of LPV and ultra LPV by reducing PaCO<sub>2</sub> and thereby normalizing pH<sup>25</sup>





## References

1. Silversides JA, Ferguson ND. *Crit Care* 2013;17:225.
2. Gattinoni L, Protti A. *CMAJ* 2008;178:1174–6.
3. Slutsky AS. *Am J Respir Crit Care Med* 2015;191:1106–15.
4. Bellani G, et al. *JAMA* 2016;315:788–800.
5. Ranieri VM, et al. *JAMA* 2012;307:2526–33.
6. Peek GJ, et al. *Lancet* 2009;374:1351–63.
7. Esteban A, et al. *Chest* 2000;117:1690–6.
8. Villar J, et al. *Crit Care Med* 2006;34:1311–8.
9. The Acute Respiratory Distress Syndrome Network. *N Engl J Med* 2000;342:1301–8.
10. Amato MB, et al. *N Engl J Med* 1998;338:347–54.
11. Amato MB, et al. *N Engl J Med* 2015;372:747–55.
12. Hess D, Kacmarek R. *Essentials of mechanical ventilation*. 3rd ed: McGraw-Hill Education; 2014.
13. Bein T, et al. *Intensive Care Med* 2013;39:847–56.
14. Fanelli V, et al. *Crit Care* 2016;20:36.
15. Dellinger RP, et al. *Crit Care Med* 2013;41:580–637.
16. World Health Organization. *Clinical Management of Severe Acute Respiratory Infections When Novel Coronavirus is Suspected*, 2013. Available at: [http://www.who.int/csr/disease/coronavirus\\_infections/InterimGuidance\\_ClinicalManagement\\_NovelCoronavirus\\_11Feb13u.pdf?ua=1](http://www.who.int/csr/disease/coronavirus_infections/InterimGuidance_ClinicalManagement_NovelCoronavirus_11Feb13u.pdf?ua=1) (Accessed July 2017).
17. Intensive Care Society. *Guidelines for the Provision of Intensive Care Services*, 2015. Available at: <http://www.ics.ac.uk/ICS/guidelines-and-standards.aspx> (Accessed July 2017).
18. Cho YJ, et al. *Tuberc Respir Dis* 2016;79:214–33.
19. Brunkhorst FM, et al. *Crit Care Med* 2008;36:2719–25.
20. Rubenfeld GD, et al. *Crit Care Med* 2004;32:1289–93.
21. Kalhan R, et al. *Crit Care Med* 2006;34:300–6.
22. Dennison CR, et al. *Crit Care Med* 2007;35:2747–54.
23. Nin N, et al. *Intensive Care Med* 2017;43:200–8.
24. Morimont P, et al. *Crit Care* 2015;19:117.
25. Terragni PP, et al. *Anesthesiology* 2009;111:826–35.
26. Del Sorbo L, et al. *Lancet Respir Med* 2014;2:154–64.