

Baxter

Acute RRT Modalities: Comparisons and Considerations





Table of Contents

1

ACUTE KIDNEY INJURY (AKI) OVERVIEW

2

RRT MODALITIES FOR AKI

3

RRT MODALITY CONSIDERATIONS

4

SUMMARY


5

ACRONYMS/ABBREVIATIONS/REFERENCES

ACUTE KIDNEY INJURY OVERVIEW



Acute kidney injury is **common** among hospitalised patients globally¹



AKI AFFECTS AN ESTIMATED **20%** OF HOSPITALISED PATIENTS WORLDWIDE^{1,*}

AKI is a **serious** condition

AKI IS ASSOCIATED WITH AN INCREASED RISK OF **MORBIDITY and MORTALITY**²⁻⁶


AKI IS ASSOCIATED WITH AN INCREASED RISK OF **CKD, including ESRD**⁷⁻⁹

*Multicentre meta-analysis of 154 studies (n=3,585,911), primarily in hospital settings, that adopted a KDIGO-equivalent AKI definition between 2004 and 2012. Pooled rates.¹


ACUTE KIDNEY INJURY OVERVIEW




Fluid overload is one condition that may **adversely** impact AKI patient prognosis^{10,11}



Consequences of fluid overload may lead to organ dysfunction¹³



Fluid overload at RRT initiation for AKI has been associated with an increased risk of mortality^{11,*}



FLUID OVERLOAD IN PATIENTS WITH **AKI** IS A **SERIOUS CONDITION**¹⁴⁻¹⁶

*Prospective, observational cohort study of 296 adults treated with RRT in 17 Finnish ICUs from Sep 2011–Feb 2012.¹¹


ACUTE KIDNEY INJURY OVERVIEW



AKI is associated with **substantial financial burden**^{17–19,*}

AKI status impacts daily costs¹⁸

PATIENTS WITH AKI HAVE SIGNIFICANTLY HIGHER DAILY COSTS COMPARED WITH PATIENTS WITHOUT AKI^{18,†}



AKI is expensive even relative to other acute medical conditions¹⁹

Acute medical condition	Adjusted mean cost difference, in 2012 USD (95% CI) ^a
AKI-D^b	11,016 (10,468, 11,564)
Sepsis	4822 (4696, 5068)
VTE	3782 (3611, 3953)
Acute pancreatitis	1802 (1676, 1929)
AKI^c	1795 (1692, 1899)
Pneumonia	1705 (1584, 1825)
Stroke	1427 (1281, 1573)
MI	14 (–91, 119)
GI bleed	–860 (–961, –759)

THE INCREMENTAL COST OF AKI-D OR AKI IS HIGHER THAN FOR MANY OTHER CONDITIONS FOUND IN HOSPITALISED PATIENTS^{19,‡}

^aCompared with reference group without the condition of interest.

^bCompared with patients without AKI. ^cIncludes patients with dialysis-requiring AKI (AKI-D).

WHILE EXPENDITURES MAY VARY BY COUNTRY,
AKI is a COSTLY CONDITION^{17–19}


*Costs for hospitalisation due to AKI may vary from country to country.

†Multicentre, retrospective cohort study of 659,945 adult hospital admissions across central China in 2013.¹⁸

‡2012 multicentre, retrospective study of 29,763,649 adult US hospitalisations without ESRD.¹⁹



Various **renal replacement modalities** are available for the **management of AKI**²⁰⁻²⁴








ACUTE RRT IS DELIVERED AS **EITHER**
A **CONTINUOUS OR**
INTERMITTENT THERAPY²⁰

RRT MODALITIES FOR AKI



Modalities **differ** in their typical characteristics²⁵

Typical RRT modality characteristics and settings for a 70-kg AKI patient^{25–27}

	CONTINUOUS THERAPIES			INTERMITTENT THERAPIES	
Parameter	CVH	CVVHD	CVVHDF	SLED*	IHD
Blood flow (Q _B , mL/min)	150–250	150–250	150–250	100–300	200–300
Predominant solute transport principle					
Ultrafiltrate (mL/h)	1500–2000	variable	1000–1500	variable	variable
Dialysate flow (Q _D , mL/h)	0	1500–2000	1000–1500	6000–18,000	18,000–30,000
Replacement fluid for zero balance (mL/h)	1500–2000	0	1000–1500	0	0
Urea clearance (mL/min)	25–33	25–33	25–33	80–90	200–500

*SLED is a type of PIRRT.²¹



Convection




Diffusion

Q_B, Q_D, AND UREA CLEARANCE TEND TO BE **LOWER** IN **CONTINUOUS** THERAPIES THAN IN **INTERMITTENT** THERAPIES^{25–27}



Individual patient needs can be addressed by considering the characteristics of the various **RRT modalities**²⁸

Relative **features, risks, and burdens** of different RRT modalities²⁸



EACH RRT MODALITY HAS POTENTIAL **BENEFITS** AND **LIMITATIONS** FOR THE MANAGEMENT OF PATIENTS WITH AKI²⁸



Selection of RRT modality requires careful consideration of many patient- and ICU-specific factors^{25,28}

Overview of modality considerations



**CLINICAL CONSIDERATIONS:
FLUID OVERLOAD AND HAEMODYNAMIC INSTABILITY**



CLINICAL CONSIDERATIONS: LONG-TERM OUTCOMES



MACHINE AND PRESCRIPTION CONSIDERATIONS



SOLUTION CONSIDERATIONS



LONG-TERM COST CONSIDERATIONS



EQUIPMENT FOOTPRINT AND MOBILITY CONSIDERATIONS

RRT MODALITY CONSIDERATIONS



Clinical considerations: fluid overload and haemodynamic instability




Fluid overload in AKI patients can be treated by fluid removal during RRT, but rapid fluid removal that does not allow time for plasma refill may lead to **haemodynamic instability**^{25,29}



Avoiding rapid fluid removal to prevent hypovolaemia may **improve** AKI patient outcomes^{25,29}

Modality comparisons³⁰



CRRT IS A **PREFERRED RRT**
BY MANY CLINICIANS FOR AKI PATIENTS WHO ARE
HAEMODYNAMICALLY UNSTABLE^{25,29}

RRT MODALITY CONSIDERATIONS




Clinical considerations: long-term outcomes



AKI is associated with an **increased risk** of long-term **dialysis dependence**;⁸ acute **RRT modality type** may impact this risk³¹⁻³⁴

Patients on chronic dialysis at day 90
by initial RRT modality^{31,*}



CHRONIC DIALYSIS HAZARD RATIO (95% CI) FOR CRRT vs IHD WAS 0.75 (0.65–0.87), $P < 0.0001$ ³¹

Modality comparisons³¹⁻³⁴

CONTINUOUS	INTERMITTENT	
CRRT	PIRRT	IHD
Patients are less likely to require chronic dialysis following initial AKI episode compared with patients treated with IHD	Insufficient evidence	It has been reported that patients are more likely to require chronic dialysis following initial AKI episode compared with patients treated with CRRT

USE OF **CRRT** FOR AKI MANAGEMENT HAS BEEN ASSOCIATED WITH A **LOWER RISK of CHRONIC DIALYSIS** COMPARED WITH IHD³¹⁻³⁴

*Retrospective multicentre cohort study of critically ill adults with AKI between 1996 and 2009. 2004 patients originally treated with CRRT and 2004 patients originally treated with IHD were propensity matched and rates of dialysis dependence were compared.³¹

RRT MODALITY CONSIDERATIONS




Machine and prescription considerations



RRT machines deliver **different** dose intensities over **different** durations of therapy^{21,22,25}

Kinetic modeling of urea clearance by different RRT modalities³⁵



A SAWTOOTH PATTERN WAS OBSERVED WHEN USING INTERMITTENT THERAPIES TO REMOVE UREA, WHILE CONTINUOUS THERAPY MAINTAINED A CONSISTENT BUN LEVEL OVER TIME³⁵

Modality comparisons^{21,22,35}

CONTINUOUS	INTERMITTENT	
CVVH	SLED	IHD
Intended to run 24 h/day	Typically run in 6–12 h sessions delivered 3–7 times/week	Typically run in 3–6 h sessions delivered 3–7 times/week
• Slow but continuous urea clearance helps avoid spikes in BUN levels	• Intermittent nature does not allow for continuous urea clearance, which could result in variable BUN levels	• Intermittent nature does not allow for continuous urea clearance, which could result in variable BUN levels

UNLIKE IHD OR PIRRT, **CRRT** IS RUN ON MACHINES THAT DELIVER **CONTINUOUS** SOLUTE REMOVAL^{22,35}

RRT MODALITY CONSIDERATIONS



Solution considerations



Typically, CRRT solutions are **commercially** prepared, while IHD and PIRRT use **local water sources** to prepare dialysate^{29,36,37}



Preparing solutions on-line from local water sources **necessitates** water **treatment** and routine water **quality monitoring** to assure clean water standards are met³⁶⁻³⁸

Modality comparisons^{29,36-40}

CONTINUOUS	INTERMITTENT	
CRRT	PIRRT	IHD
<p>Because no on-line solutions are typically used, no water treatment systems are required</p> <ul style="list-style-type: none">• <p>Monitoring water quality is not applicable</p>	<p>If a centralized water treatment system is unavailable in the ICU, individual water quality monitoring is necessary</p> <ul style="list-style-type: none">• <p>If a centralized water treatment system is not used, staff need to monitor dialysate quality for individual patients</p> <ul style="list-style-type: none">• <p>Disinfection requirements may limit treatment duration to <12 hours⁴¹</p>	<p>If a centralized water treatment system is unavailable in the ICU, individual water quality monitoring is necessary</p> <ul style="list-style-type: none">• <p>If a centralized water treatment system is not used, staff need to monitor dialysate quality for individual patients</p>

WATER TREATMENT AND QUALITY TESTING MAY CONTRIBUTE
TO **INCREASED MONITORING**
WHEN USING SOLUTIONS PREPARED ON-LINE FOR
IHD and PIRRT^{39,42}

RRT MODALITY CONSIDERATIONS




Long-term cost considerations



Because **initial RRT modality** may impact the risk of **chronic dialysis**,³¹ **long-term costs** of AKI may also be influenced by **initial treatment modality**⁴³

Cumulative costs of dialysis dependence by initial AKI treatment modality^{43,*}



MEAN 5-YEAR TOTAL COST/PATIENT OF AKI-D⁺ WAS \$37,780 FOR CRRT AS THE INITIAL MODALITY COMPARED WITH \$39,448 FOR IRRT⁴³

[†]Including cost of dialysis dependence. Cost in 2013 USD.

Modality comparisons⁴³

CONTINUOUS		INTERMITTENT	
CRRT		PIRRT	IHD
Total costs may be lower due, in part, to a lower risk of chronic dialysis		Insufficient evidence to compare to CRRT or IHD	Total costs may be higher due, in part, to a higher risk of chronic dialysis

THE **LONG-TERM COST** OF AKI MAY BE **LOWER** FOR PATIENTS INITIALLY TREATED WITH **CRRT** COMPARED TO THOSE TREATED WITH IHD⁴³

*Health outcomes and healthcare costs were simulated and averaged for a cohort of 1000 patients initiated on CRRT and a cohort of 1000 patients initiated on IRRT. All costs were inflated to 2013 USD.⁴³



Equipment footprint and mobility considerations



Water treatment systems required for IHD and PIRRT add to physical **space** requirements and water lines may limit RRT mobility in ICUs without central water treatment systems^{37,40}



In ICUs without central water treatment systems, portable **water treatment devices** may be necessary,⁴⁰ which can occupy as much as 0.13–0.16 m² of floor space^{44,45}

Modality comparisons^{36,37,40,41,46–49}

CONTINUOUS	INTERMITTENT	
	PIRRT	IHD
CRRT	<p>Both the IHD machine and water treatment systems contribute to the therapy's physical footprint, which may impact treatment mobility in ICUs without central water treatment systems</p> <ul style="list-style-type: none"> In situations where a central water treatment system is not utilised, the greater physical footprint of the machine + water treatment system may impact ICU spacing 	

Because the CRRT machine is the **only component** that contributes to the therapy's physical footprint, treatment mobility may be increased

- No **space considerations** for water treatment systems are necessary

WATER TREATMENT EQUIPMENT MAY ADD TO THE FOOTPRINT OF **IHD** AND **PIRRT** SYSTEMS, POTENTIALLY DECREASING TREATMENT MOBILITY AND IMPACTING SPACING CONSIDERATIONS^{40,47–49}

SUMMARY



AKI is a **common** and **costly** condition among ICU patients,^{1,17-19} and is associated with increased risks of **morbidity and mortality**²⁻⁹



Acute RRT is delivered as **either a continuous or an intermittent** therapy, each of which have unique characteristics, settings, and limitations^{20,25-28}

Selection of RRT modality requires careful consideration of many patient- and ICU-specific factors^{25,28}



FLUID OVERLOAD AND
HAEMODYNAMIC
INSTABILITY



LONG-TERM
CLINICAL
OUTCOMES



MACHINE AND
PRESCRIPTION



SOLUTIONS



LONG-TERM COSTS



EQUIPMENT
FOOTPRINT AND
MOBILITY

CRRT IS A PREFERRED
RENAL REPLACEMENT THERAPY
BY MANY CLINICIANS FOR PATIENTS WITH AKI
WHO ARE HAEMODYNAMICALLY UNSTABLE^{25,28}





AKI, acute kidney injury; AKI-D, dialysis-requiring AKI; BUN, blood urea nitrogen; CI, confidence interval; CKD, chronic kidney disease; CRRT, continuous renal replacement therapy; CVVH, continuous veno-venous haemofiltration; CVVHD, continuous veno-venous haemodialysis; CVVHDF, continuous veno-venous haemodiafiltration; dL, decilitre; EDHF, extended daily haemofiltration; ESRD, end-stage renal disease; Feb, February; GI, gastrointestinal; h, hour; HD, haemodialysis; ICU, intensive care unit; IHD, intermittent haemodialysis; IRRT, intermittent renal replacement therapy; K⁺, potassium ion; KDIGO, Kidney Disease Improving Global Outcomes; kg, kilogram; m², square meters; mg, milligram; MI, myocardial infarction; min, minute; mL, millilitre; PD, peritoneal dialysis; PIRRT, prolonged intermittent renal replacement therapy; NH₃, ammonia; Q_B, blood flow rate; Q_D, dialysis flow rate; RRT, renal replacement therapy; Sep, September; SLED, sustained or slow low-efficiency dialysis; SLEDD, sustained or slow low-efficiency daily dialysis; US, United States; USD, United States dollar; vs, versus; VTE, venous thromboembolism

1. Susantitaphong P, et al. *Clin J Am Soc Nephrol*. 2013;8:1482-1493.
2. Coca SG, et al. *Am J Kidney Dis*. 2009;53:961-973.
3. Ricci Z, et al. *Kidney Int*. 2008;73:538-546.
4. Chawla LS, et al. *Clin J Am Soc Nephrol*. 2014;9:448-456.
5. Brown JR, et al. *Ann Thorac Surg*. 2016;102:1482-1489.
6. Wu VC, et al. *J Am Soc Nephrol*. 2014;25:595-605.
7. Ishani A, et al. *J Am Soc Nephrol*. 2009;20:223-228.
8. Coca SG, et al. *Kidney Int*. 2012;81:442-448.
9. Wald R, et al. *JAMA*. 2009;302:1179-1185.
10. Bellomo R, et al. *Crit Care Med*. 2012;40:1753-1760.
11. Vaara ST, et al. *Crit Care*. 2012;16:R197.
12. Kaddourah A, et al. *N Engl J Med*. 2017;376:11-20.
13. Prowle JR, et al. *Nat Rev Nephrol*. 2010;6:107-115.
14. Zhang L, et al. *J Crit Care*. 2015;30:860.e7-13.
15. Bouchard J, et al. *Kidney Int*. 2009;76:422-427.
16. Heung M, et al. *Nephrol Dial Transplant*. 2012;27:956-961.
17. National Clinical Guideline Centre (UK). Acute Kidney Injury: Prevention, Detection and Management Up to the Point of Renal Replacement Therapy. 2013; NICE Clinical Guidelines, No. 169. Introduction.
18. Xu X, et al. *Clin J Am Soc Nephrol*. 2015;10:1510-1518.
19. Silver SA, et al. *J Hosp Med*. 2017;12:70-76.
20. Fleming GM. *Organogenesis*. 2011;7:2-12.
21. O'Reilly P, Tolwani A. *Crit Care Clin*. 2005;367-378.
22. Pannu N, Gibney RTN. *Ther Clin Risk Manag*. 2005;1:141-150.
23. Sun Z, et al. *Crit Care*. 2014;18:R70.
24. Kitchlu A, et al. *BMC Nephrol*. 2015;16:127.
25. Kidney Disease: Improving Global Outcomes (KDIGO) Acute Kidney Injury Work Group. KDIGO Clinical Practice Guideline for Acute Kidney Injury. *Kidney Int Suppl*. 2012;2:1-138.
26. Baxter Corporation. REVACLEAR Dialyzer Technology. 2017. Available from: http://www.baxter.ca/en_CA/assets/downloads/2017/Revaclear%20Spec%20Sheet%20Brochure%20English.pdf (accessed December 2018).
27. Fresenius Medical Care. Optiflux Dialyzers. 2016. Available from: http://www.fmcna-dialyzers.com/images/pdf/101046-Optiflux-Dialyzer_SpecSheet.pdf (accessed December 2018).
28. Ostermann M, et al. *Blood Purif*. 2016;42:224-237.
29. Murugan R, et al. *Blood Purif*. 2016;42:266-278.
30. Ostermann M. In: Bellomo R, et al (eds). 40 Years of Continuous Renal Replacement Therapy. *Contrib Nephrol*. 2018;194:51-59.
31. Wald R, et al. *Crit Care Med*. 2014;42:868-877.
32. Bell M, et al. *Intensive Care Med*. 2007;33:773-780.
33. Cartin-Ceba R, et al. *Intensive Care Med*. 2009;35:2087-2095.
34. Lin YF, et al. *Am J Surg*. 2009;198:325-332.
35. Liao Z, et al. *Artif Organs*. 2003;27:802-807.
36. Coulliette AD, Arduino MJ. *Semin Dial*. 2013;26:427-438.
37. Glorieux G, et al. *Nephrol Dial Transplant*. 2012;27:4010-4021.
38. Azar AT, Ahmad S. Hemodialysis Water Treatment System. In: Azar A (eds). *Modelling and Control of Dialysis Systems. Studies in Computational Intelligence*. 2013;404:347-378.
39. BC Renal Agency. Clinical Practice Standards and Procedures for Dialysis Water Quality: 2b: Endotoxin Testing of Dialysis Water. 2011. Available from: http://www.bcrenalagency.ca/resource-gallery/Documents/2bEndotoxin-Testing-of-Dialysis-Water-Final_2012.pdf (accessed December 2018).
40. Bellomo R, et al. *Crit Care Resusc*. 2002;4:281-290.
41. Fresenius Medical Care. 2008K² Hemodialysis Machine Operator's Manual. 2016. Available from: https://fmcna.com/wp-content/uploads/documents/490136_Rev_J.pdf (accessed December 2018).
42. Food and Drug Administration. Quality Assurance Guidelines for Hemodialysis Devices. 1991. Available from: <https://www.fda.gov/downloads/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm073435.pdf> (accessed December 2018).
43. Ethgen O, et al. *Nephrol Dial Transplant*. 2015;90:54-61.
44. Mar Cor Purification. 700 Series Portable. 2014. Available from: [http://www.mcpur.com/main/library/12_brochures/1238642_\(700\).pdf](http://www.mcpur.com/main/library/12_brochures/1238642_(700).pdf) (accessed December 2018).
45. Mar Cor Purification. Millenium HX. 2015. Available from: [http://www.mcpur.com/main/library/12_brochures/3027573_\(MHX\).pdf](http://www.mcpur.com/main/library/12_brochures/3027573_(MHX).pdf) (accessed December 2018).
46. Baxter Healthcare Corporation. The PRISMAFLEX System. 2016. Available from: https://www.baxter.com/sites/g/files/ebysai746/files/2017-11/Prismaflex-07.11-Brochure-New_Accts.pdf (accessed December 2018).
47. Ledebro I, Blankestijn PJ. *NDT Plus*. 2010;3:8-16.
48. Kasperek T, Rodriguez OE. *Clin J Am Soc Nephrol*. 2015;10:1061-1071.
49. Poeppel K, et al. *Vet Clin Small Anim*. 2011;177-191.