

Polyflux L

DESIGNED FOR:

LFHD [Low flux]

MEMBRANE:

POLYAMIX [PAES/PVP/PA, BPA-free]

THE PROVEN BALANCE OF QUALITY AND PERFORMANCE IN LOW-FLUX

The **Polyflux L** dialyzer series is specialized for low-flux hemodialysis treatments, featuring a distinctive membrane acting as an effective barrier to potential fluid contaminants¹, while still delivering high performance². **Polyflux L** dialyzers are a good choice for proven biocompatible yet effective low-flux therapies, designed with safety in mind.

DESIGNED TO PROMOTE BIOCOMPATIBILITY²

The **Polyflux L** dialyzers are designed to deliver high-quality low-flux hemodialysis treatments.

- Since 1988, over 300 million **Polyflux** dialyzers have been used globally³
- The **Polyflux L** dialyzers are designed to prevent endotoxins from crossing the dialyzer membrane^{1,2}
- The **Polyflux L** dialyzers are steam sterilized inside-out, to promote biocompatibility, avoiding exposure to chemicals such as ethylene oxide and manufacturing residues^{4,5}

WITH HIGH PERFORMANCE IN MIND

The **Polyflux L** dialyzers feature an exclusive 3-layered membrane structure, designed to support a stable high performance over time.

- Effective clearance of standard dialysis markers, such as urea or phosphates⁶
- A clinical case series study suggests the **Polyflux L** dialyzers may reduce the signs and symptoms of hemodialysis-associated eosinophilia⁷



Polyflux L Specifications

MATERIALS	POLYFLUX 14 L	POLYFLUX 17 L	POLYFLUX 21 L
Membrane	Polyamix Polyarylethersulfone, Polyvinylpyrrolidone and Polyamide blend BPA-free		
Potting	Polyurethane (PUR)		
Housing	Polycarbonate (PC)		
Gaskets	Silicone rubber (SIR)		
Protection caps	Polypropylene (PP)		
Sterilization	Steam (inside-out)		
Sterile barrier	Medical Grade Paper		
SPECIFICATIONS			
UF-Coefficient (mL/h*mmHg)*	10	12.5	15
KoA urea*	851	1026	1268
Blood Compartment volume (mL)	81	104	123
Minimum recommended priming volume (mL)	500		
Maximum TMP (mmHg)	600		
Recommended Q _B (mL/min)	200-400	200-500	300-500
Storage conditions	<30°C (or <86°F)		
Units per box	24		
Gross/net weight (g)	254/225	274/245	294/265
MEMBRANE			
Effective Membrane Area (m ²)	1.4	1.7	2.1
Fiber inner diameter (µm)	215		
Fiber wall thickness (µm)	50		

* According to EN 1283/ISO 8637:

- UF-Coefficient: measured with bovine blood, Hct 32%, Pct 60g/L, 37°C
- KoA urea: calculated at Q_B=300 mL/min, Q_D=500mL/min, UF=0 mL/min
- Clearances In-Vitro: measured at UF=0 mL/min, ±10%

CLEARANCES IN VITRO (mL/min)*	POLYFLUX 14 L	POLYFLUX 17 L	POLYFLUX 21 L
Urea (60 Da) (Q_B-Q_D, mL/min)			
200/500	190	194	
300/500	252	264	275
400/500	293	310	328
500/500		342	364
200/700	194	197	
300/700	267	276	285
400/700	319	336	343
500/700		380	403
Phosphate (95 Da)			
200/500	152	163	
300/500	183	200	218
400/500	203	224	247
500/500		240	267
200/700	160	170	
300/700	197	213	231
400/700	221	242	266
500/700		264	272
Creatinine (113 Da)			
200/500	171	179	
300/500	214	230	246
400/500	241	262	283
500/500		284	310
200/700	178	185	
300/700	229	244	258
400/700	264	284	306
500/700		313	341
Vitamin B12 (1.4 kDa)			
200/500	90	101	
300/500	100	114	131
400/500	106	122	142
500/500		128	149
200/700	96	107	
300/700	107	121	138
400/700	114	130	150
500/700		137	159

1. Schepers E, Glorieux G, Elout S, et al. *Assessment of the association between increasing membrane pore size and endotoxin permeability using a novel experimental dialysis simulation set up.* BMC Nephrology. 2018; 19:1.
2. Ronco C, et al. *Evolution of synthetic membranes for blood purification: the case of the Polyflux family.* Nephrol Dial Transplant 2003;18(Suppl 7):viii10-20.
3. Baxter. Data on file. *Dialyzers Sales Report.* 2018.
4. Golli-Bennour EE, et al. *Cytotoxic effects exerted by polyarylsulfone dialyzer membranes depend on different sterilization processes.* Int Urol Nephrol. 2011; 43:483-490.
5. D'Ambrosio FP, et al. *Ethylene oxide allergy in dialysis patients.* Nephrol Dial 1997;12:1461-1463.
6. Krause B, et al. *Polymeric Membranes for Medical Applications.* Chemie Ingenieur Technik. 2003; 75:1725-1732.
7. Tielemans C, et al. *Clinical assessment of Performance and Blood Compatibility Profile of a New Synthetic Low Flux Hemodialyzer.* Blood Purif. 2002; 20:214-215.

The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (Full quality assurance system of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/47/EC)

For safe and proper use of the device, please refer to the Instructions for Use

CE 0086

Baxter, Polyamix and Polyflux are trademarks of Baxter International Inc. or its subsidiaries.

GLBL/MG3/18-0020f - January 2019

MANUFACTURER
Gambro Dialysatoren GmbH
Holger-Crafoord-Strasse 26
72379 Hechingen
Germany

Baxter Healthcare Corporation
One Baxter Parkway
Deerfield, IL 60015
USA
1-800-422-9837