

**Baxter**

Novaline

BLOODLINE TUBING SYSTEM

# BLOODLINE TUBING SYSTEM

# SIMPLE

A clear, flexible medical tubing system is shown against a light blue background with geometric patterns. The tubing is coiled into a large loop. A red stopcock is attached to the tubing in the upper part of the loop. At the bottom of the loop, there is a red connector. The tubing ends in a clear, flexible connector with a red ring. The overall design is clean and professional.

# EFFICIENT

# MADE FOR YOU AND YOUR PATIENTS

Over the years, we have worked in close cooperation with dialysis clinics, to develop product features that help meeting the requirements of dialysis providers and patients.



## EASY HANDLING AND TRACEABILITY

A clinician-requested feature of our bloodlines is the removable packaging label, which can be peeled off and attached to the patient records.

## COMBI-CONNECTOR

Allows for an easy placement of the blood pump segment into the blood pump.




## DUAL MATERIAL DIALYZER CONNECTOR

A dialyzer connector design that supports a robust connection to the dialyzer, to help minimize the risk of leaks.



## PRESSURE TRANSDUCER PROTECTOR

There is a single or double pressure transducer protector housing, which contains one membrane made of two layers.



## NOVALINE MODELS FOR HD ON AK SYSTEMS

CODE	MODEL	SYSTEM	DN/SN	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE-DIALYZER CHAMBER	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955299	BL 05	AK 95 S, AK 96, AK 98	DN	24	Eto	DEHP	8x12 mm 250 mm	Yes	No	1	Yes	1	Mono	No
955419	BL 05R	AK 95 S, AK 96, AK 98	DN	24	Ebeam	DEHP Free	8x12 mm 250 mm	Yes	No	1	Yes	1	Dual	No
955491	BL 06	AK 95 S, AK 96, AK 98	DN	24	Eto	DEHP	8x12 mm 250 mm	Yes	No	1	Yes	1	Mono	No
955301	BL 10	AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	24	Eto	DEHP	8x12 mm 250 mm	No	No	1	No	2	Mono	No
955300	BL 10 R	AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	24	Radiation	DEHP Free	8x12 mm 250 mm	No	No	1	No	2	Dual	No
955304	BL 14	AK 95 S, AK 96, AK 98, AK 200 S	DN	24	Eto	DEHP	8x12 mm 250 mm	Yes	No	1	No	2	Dual	No
955307	BL 24	AK 95 S, AK 96, AK 98, AK 200 S	DN	24	Radiation	DEHP Free	8x12 mm 280 mm	Yes (pre-attached)	No	1	Yes	1	Dual	Yes
955308	BL 25	AK 95 S, AK 96, AK 98, AK 200 S	DN	24	Radiation	DEHP Free	8x12 mm 280 mm	Yes (pre-attached)	No	1	No	1	Mono	No
955312	BL 90 D4	AK 95 S, AK 96, AK 98	DN	24	Ebeam	DEHP Free	8x12 mm 250 mm	No	No	1	Yes	1	Dual	No
955462	BL 95 SN	AK 95 S, AK 96, AK 98	SN SP	12	Radiation	DEHP Free	8x12 mm 250 mm	Yes	No	1	Yes	2	Dual	Yes
955798	BL 102*	AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	24	Radiation	DEHP free	8x12 mm 250 mm	Yes	No	1	No	1	Dual	Yes
955303	BL 105	AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	20	Radiation	DEHP Free	8x12 mm 250 mm	No	No	1	Yes	2	Dual	Yes
955443	BL 106	AK 95 S, AK 96, AK 98,	DN	20	Radiation	DEHP Free	8x12 mm 250 mm	No	No	2	Yes	2	Dual	Yes
955305	BL 148 SN	AK 200 S, AK 200 Ultra S	SN DP	12	Radiation	DEHP Free	P1: 8x12 mm 250 mm P2: 8x12 mm 250 mm	No	No	1	Yes	2	Dual	No
955448	BL 207	AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	24	Radiation	DEHP Free	8x12 mm 247 mm	No	Yes	1	No	2	Dual	Yes
955445	BL 208	AK 200 S, AK 200 Ultra S	DN	20	Radiation	DEHP Free	8x12 mm 247 mm	No	Yes	2	Yes	2	Dual	Yes
955447	BL 211 SN	AK 200 S, AK 200 Ultra S	SN DP	12	Radiation	DEHP Free	P1: 8x12 mm 247 mm P2: 8x12 mm 250 mm	No	Yes	2	Yes	2	Dual	Yes
955446	BL 224	AK 200 S, AK 200 Ultra S	DN	20	Radiation	DEHP Free	8x12 mm 247 mm	Yes (pre-attached)	Yes	2	Yes	2	Dual	Yes
955449	BL 245	AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S	DN	24	Radiation	DEHP free	8x12 mm 247 mm	Yes	Yes	1	No	1	Dual	Yes

\* Only for Canada

## NOVALINE MODELS FOR HD INTEGRA SYSTEM

CODE	MODEL	SYSTEM	DN/SN	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE-DIALYZER CHAMBER	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955559	HINT-4DG	Integra/ Integra Alpha	DN	20	Gamma	DEHP Free	8x12 mm 340 mm	No	No	2	Yes	2	Dual	Yes
955560	HINT-11G	Integra	SN DN	12	Gamma	DEHP Free	P1: 8x12 mm 340 mm P2: 8x12 mm 340 mm	No	No	2	Yes	2	Dual	Yes

## NOVALINE MODELS FOR HD, HF AND HDF ON AK SYSTEMS

CODE	MODEL	SYSTEM	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE-DIALYZER CHAMBER	PRE-DILUTION	POST-DILUTION	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955302	BL 10 E2C	AK 200 S, AK 200 Ultra S	24	Radiation	DEHP Free	8x12 mm 250 mm	No	No	1	Yes	Yes	Yes	2	Dual	No
955313	BL 90 D4T	AK 200 S, AK 200 Ultra S	24	Ebeam	DEHP Free	8x12 mm 260 mm	No	No	1	Yes	Yes	Yes	2	Dual	No
955306	BL 150 D	AK 200 S, AK 200 Ultra S	20	Radiation	DEHP Free	8x12 mm 250 mm	No	No	1	Yes	Yes	Yes	2	Dual	No
955444	BL200 HDF	AK 200 S, AK 200 Ultra S	20	Radiation	DEHP Free	8x12 mm 247 mm	No	Yes	2	Yes	Yes	Yes	2	Dual	Yes

## NOVALINE FOR HD ON OTHER MONITORS

CODE	MODEL	SYSTEM	DN/SN	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE-DIALYZER CHAMBER	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
95541	BL 31	Braun DIALOG +	DN	24	Eto	DEHP	6.3x9.8 mm 370mm	No	No	1	Yes	1	Dual	No
955451	BL 31 R	Braun DIALOG +	DN	24	Radiation	DEHP Free	8x12 mm 374 mm	Yes	No	1	Yes	1	Dual	No
955664	BL 32 R	Fresenius 4008	DN	24	Radiation	DEHP Free	8x12 mm 350 mm	No	No	1	Yes	2	Dual	No
955440	BL 33	Fresenius 4008	DN	24	Radiation	DEHP Free	8x12 mm 400 mm	Yes	No	1	Yes	1	Dual	No
955310	BL 36	Fresenius 4008	DN	24	Eto	DEHP	8x12 mm 350 mm	No	No	1	Yes	1	Mono	No
955309	BL 36 R	Fresenius 4008	DN	24	Radiation	DEHP Free	8x12 mm 350 mm	No	No	1	Yes	1	Mono	No
955459	BL 37	Fresenius 4008	DN	20	Radiation	DEHP Free	8x12 mm 355 mm	Yes	No	1	Yes	2	Dual	No
955460	BL 38 SN	Fresenius 4008	SN DP	16	Radiation	DEHP Free	P1: 8x12 mm 355 mm P2: 8x12 mm 355 mm	No	No	1	Yes	1	Dual	No
955441	BL 39	Nikkiso DBB-05-07	DN	20	Radiation	DEHP Free	8x12 mm 260 mm	Yes (pre-attached)	No	1	Yes	2	Dual	No
955442	BL 40 SN	Nikkiso DBB-05-07	SN DP	12	Radiation	DEHP Free	P1: 8x12 mm 260 mm P2: 8x12 mm 260 mm	No	No	1	Yes	1	Dual	No
955450	BL 41	Nikkiso DBB-06	DN	24	Radiation	DEHP Free	8x12 mm 265 mm	Yes (pre-attached)	No	1	Yes	1	Mono	No

## NOVALINE – LOW VOLUME AND VERY LOW VOLUME ON AK SYSTEMS

CODE	MODEL	SYSTEM	DN/SN	QUANTITY PER BOX	TOTAL VOLUME CC	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE-DIALYZER CHAMBER	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955463	BL 120N	<b>AK 98, AK 200 S, AK 200 Ultra S</b>	DN	36	32	Radiation	DEHP free	3.9 x 7.9 mm 250 mm	No	No	1	No	2	Dual	Yes
955774	BL 120N*	<b>AK 98, AK 200 S, AK 200 Ultra S</b>	DN	36	32	Radiation	DEHP free	3.9 x 7.9 mm 250 mm	No	No	1	No	2	Dual	Yes
955464	BL 121P HDF	<b>AK 98, AK 200 S, AK 200 Ultra S</b>	DN	24	87	Radiation	DEHP free	6.35x10.35 mm 247 mm	No	Yes	1	No	2	Dual	Yes
955457	BL 122P	<b>AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S</b>	DN	24	116	Radiation	DEHP free	6.35x10.35 mm 250 mm	No	No	1	Yes	2	Dual	Yes
955775	BL 122P*	<b>AK 95 S, AK 96, AK 98, AK 200 S, AK 200 Ultra S</b>	DN	24	116	Radiation	DEHP free	6.35x10.35 mm 250 mm	No	No	1	Yes	2	Dual	Yes
955461	BL 123P SN	<b>AK 200 S, AK 200 Ultra S</b>	SN DP	12	157	Radiation	DEHP free	P1: 6.35x10.35 250 mm P2: 6.35x10.35 250 mm	No	No	1	Yes	1	Dual	Yes


\* Only for Canada

DN: Double needle  
SN: Single needle  
DP: Double pump  
Radiation = gamma or beta

Please check with your local sales office for the latest information about the bloodline models available.

## ACCESSORIES


### RE-INFUSION LINE



**I202**  
105890

**Sterilization: EtO  
DEHP-free**

Infusion line with unvented spike (200 cm, with Male/Female connection) to use with close patient connection, between the bloodline arterial line and the patient access



**I201**  
105889

**Sterilization: EtO  
DEHP-free**


Infusion line with air vented spike (200 cm, with Male/Female connection) to use with close patient connection, between the bloodline arterial line and the patient access



**SP-339G**  
107472

**Sterilization: E-beam  
DEHP-free**


Infusion line (100 cm, with Male connection) for priming/saline infusion



**SP-235**  
103223

**Sterilization: EtO  
DEHP-free**

Infusion line (120 cm, with Male/Male connection) for priming/saline infusion



**MH277/SP332G**  
115568

For the **Integra** system

**Sterilization: E-beam  
DEHP-free**


Re-infusion line for AFB treatment/ low volume HDF with the **Integra** system, with spikes

NOTE: all the accessories are LATEX FREE



# ACCESSORIES


## CONNECTOR



**EVACLEAN CONNECTOR**  
114376

Sterilization: EtO

It may be used either as a recirculation adapter during use to short-circuit the arterial and venous patient lines or as a means to connect the blood tubing set to the Waste Handling Option (WHO®) on Innova and Artis/Evosys machines during the priming phase



**S-660-C**  
101354

Sterilization: EtO  
DEHP-free

Y-connector (10 cm lines, with Male/2 Female connections) to use with for example single lumen catheter




**S-1016**  
101023

Sterilization: EtO  
DEHP-free

Extension line Male/Female luer (20 cm) to extend the length of the bloodline patient line

## SAMPLING BAG/KIT




**GMB-SP06**  
115648

For Artis system, Artis Physio system, the AK 200 system, and AK 200 Ultra system

Sterilization: EtO  
DEHP-free

Sampling bag 500 ml for micro-biological or chemical dialysis fluid sampling with Male connection



**SP-427G**  
113226

For the AK 200 system, and AK 200 Ultra system

Sterilization: E-beam

Sampling kit of dialysate bath to check its microbiologic quality

## ADAPTOR



**SP-27**  
101957

Sterilization: EtO  
DEHP-free

Connector Female/Female luer (8 cm) to be used for example when pause treatment (recirculation) is required

## PRESSURE TRANSDUCER



**SP-221**  
100583

Sterilization: EtO  
DEHP-free

Line with pressure transducer protector (75 cm, with Male/Female connection) when a transducer protector needs to be changed/added



**SP-C35**  
106024

Sterilization: EtO  
DEHP-free

Pressure transducer protector (rigid connection, with Male/Female connection) when only the protector needs to be changed/added


## SPIKE



**SP-176**  
101422

Sterilization: EtO  
DEHP-free

Spike air vented when bottles are used




**SP-C13**  
101593

Sterilization: EtO  
DEHP-free

Spike unvented with Female connection when bags are used

# ACCESSORIES

## DRAIN/RINSING BAG




**SP-C36**  
100922

Sterilization: EtO  
DEHP-free

3L priming waste bag

## HEPARIN LINE




**H102C**  
105767

Sterilization: EtO  
DEHP-free

Heparin line (110 cm, with  
Male/Female connection),  
regional heparinization

## EXPANSION CHAMBER



**C705**

101739  
For the **AK 200** system,  
and **AK 200 Ultra** system

Sterilization: EtO  
DEHP-free

Expansion chamber for  
emergency switch from double  
needle to single needle/single  
pump mode

NOTE: all the accessories are LATEX FREE

For safe and proper use of products mentioned herein, please refer to the appropriate Instructions for Use.

Accessories

Novaline

CE 0086

CE 0197

Baxter, AK 95 S, AK 96, AK 200 S, AK 200 ULTRA S, Artis,  
Artis Physio, Integra and Novaline are trademarks of  
Baxter International Inc. or its subsidiaries.

4008 is a trademark of Fresenius Medical Care Deutschland GmbH.

DBB is a trademark of Nikkiso.

Dialog is a trademark of B Braun Avitum AG.

GLBL/G69/18-0034a(1) – July 2019

[renalcare.baxter.com](http://renalcare.baxter.com)  
Baxter Healthcare Corporation  
One Baxter Parkway  
Deerfield, IL 60015  
USA  
1-800-422-9837