

# BICART SELECT SYSTEM



- Individualizing dialysis fluid composition to meet the needs of the patient
- Ensures delivery of electrolytes and buffer as prescribed
- The lightweight design allows for easy handling, storage and logistics



# BICART SELECT SYSTEM SPECIFICATIONS

SPECIFICATIONS	SELECTBAG ONE CONTAINER	BICART SELECT COMBI-PAK
Intended use	<b>SelectBag One</b> is intended to be used as an acetic acid based concentrate together with the <b>BiCart</b> cartridge and the <b>SelectCart</b> cartridge for on-line preparation of hemodialysis, hemodiafiltration and hemofiltration fluids on Gambro dialysis monitors.	The <b>BiCart</b> and <b>SelectCart</b> products are intended for preparation of hemodialysis solutions and, when used on compatible Gambro dialysis systems, hemodiafiltration and hemofiltration substitution fluids. The <b>BiCart</b> and <b>SelectCart</b> products must always be used together with a suitable A-component concentrate.
Content	1 L of acid-based concentrate to be diluted	<b>BiCart</b> cartridge: 720g of NaHCO <sub>3</sub> <b>SelectCart</b> cartridge: 1200 g NaCl
Dilution Ratio	1:200	
Volume of Dialysis fluids when mixed	200 L	190 L at Na <sup>+</sup> /HCO <sub>3</sub> setting 140/34
Duration time at Q <sub>d</sub> of 500 ml/min	6H40	<b>BiCart</b> cartridge: 6H45 <b>SelectCart</b> cartridge: 6H15
Container Material	Low density polyethylene (LDPE)	Polypropylene (PP)
Packaging	Cardboard box with 6 bags	Cardboard box with 8 cartridges
Box size	H 320 mm x L 258 mm x W 182 mm	H 298 mm x L 198 mm x W 322 mm
Box weight	Approximately 7 kg	Approximately 9 kg
Quantity per pallet	300 bags	384 Cartridges
Shelf life	24 months	24 months
Storage condition	Store between +4°C and +30°C	Store below +40°C

SelectBag One		Contribution from A to the dialysis fluid				
PRODUCT CODE	PRODUCT NAME	K+ mmol/l	CA <sup>2+</sup> mmol/l	MG <sup>2+</sup> mmol/l	ACETATE mmol/l	GLUCOSE g/l
112629	AX125G	1.0	1.25	0.5	3	1
112630	AX150G	1.0	1.50	0.5	3	1
112631	AX175G	1.0	1.75	0.5	3	1
112634	AX225G	2.0	1.25	0.5	3	1
110532	AX250G	2.0	1.50	0.5	3	1
112635	AX260G	2.0	1.60	0.5	3	1
112636	AX275G	2.0	1.75	0.5	3	1
112641	AX325G	3.0	1.25	0.5	3	1
112642	AX350G	3.0	1.50	0.5	3	1
112644	AX375G	3.0	1.75	0.5	3	1
112646	AX425G	4.0	1.25	0.5	3	1
112647	AX450G	4.0	1.50	0.5	3	1
112648	AX475G	4.0	1.75	0.5	3	1

SelectBag One		Contribution from A to the dialysis fluid				
PRODUCT CODE	PRODUCT NAME	K+ mmol/l	CA <sup>2+</sup> mmol/l	MG <sup>2+</sup> mmol/l	ACETATE mmol/l	GLUCOSE g/l
112621	AX125	1.0	1.25	0.5	3	0
112622	AX150	1.0	1.50	0.5	3	0
112623	AX225	2.0	1.25	0.5	3	0
110531	AX250	2.0	1.50	0.5	3	0
112624	AX275	2.0	1.75	0.5	3	0
112625	AX325	3.0	1.25	0.5	3	0
112626	AX350	3.0	1.50	0.5	3	0

BiCart Select combi-pak		Contribution from A to the dialysis fluid		
PRODUCT CODE		SODIUM mmol/l	CHLORIDE mmol/l	BICARBONATE mmol/l
955847 EEA / 955848 WW / 955852 CHINA		140	107-111	34

For safe and proper use of the device, please refer to the Instructions for Use

1. Package leaflet.

